

## A CLEANER & DISINFECTANT PRODUCT FOR EVAPORATOR COILS IN THE HVACR INDUSTRY

- Cleans & Deodorizes in one easy step
- EPA registered for the HVAC application
- Ready-to-use disinfectant formulation
- No rinsing required
- Eliminates odors
- Cuts dirt and grime

### Description

Evap-Fresh is a ready-to-use, non-abrasive, one-step cleaner and disinfectant for evaporator coils in small air conditioners, commercial air handling equipment, HVAC cooling coils and refrigeration equipment. It is the first product of its kind in the industry to receive EPA registration for such claims.

Evap-Fresh has been tested effective via AOAC methodology against a wide variety of organisms:

- Salmonella typhi
- Yersinia enterocolitica
- Enterococcus faecium
- Corynebacterium ammoniagenes
- Methicillin resistant Staphylococcus aureus (MRSA)
- Vancomycin resistant Enterococcus faecalis
- Methicillin resistant Staphylococcus epidermidis
- Vancomycin intermediate resistant Staphylococcus aureus
- Pseudomonas aeruginosa
- Staphylococcus aureus
- Salmonella choleraesuis
- Escherichia coli
- Streptococcus pyogene
- Listeria monocytogenes

## Indoor Air Quality Products

### Evap Fresh™



### Directions for Use

#### Gallon:


Hold container six to eight inches from surface to be treated. Spray area until it is covered with the solution. Allow product to penetrate and remain wet for **3 minutes**. No scrubbing necessary. The product will not leave grit or soap scum.

#### Aerosol:

Shake can before using. Hold container six to eight inches from surface to be treated. Spray area until it is covered with white foam. Allow product to penetrate and remain wet for **10 minutes**. No scrubbing necessary. The product will not leave grit or soap scum.

### Packaging

- |                      |                |
|----------------------|----------------|
| 1 gallon             | <b>4166-08</b> |
| 18 oz. aerosol spray | <b>4166-75</b> |



## Addition Information

- Exhibits effective disinfectant activity against *Clostridium difficile* vegetative form, AATC 9689, when the treated surface is allowed to remain wet for 1 minute (60 seconds).
- Tuberculocidal Activity - Evap-Fresh exhibits disinfectant efficacy against *Mycobacterium tuberculosis* (BCG) at 20 degrees centigrade when the treated surface is allowed to remain wet for 5 minutes.
- Virucidal Activity - Evap-Fresh kills on hard, non-porous inanimate surfaces when allowed to remain wet for a 10 minute contact time against: Hepatitis A Virus (HAV), Canine Parvovirus and Poliovirus Type 1.
- When used on environmental, inanimate hard surfaces exhibits virucidal activity against Hepatitis B Virus (HBV), Hepatitis C Virus (HCV), and Bovine Viral Diarrhea Virus (BVDV) when the surface is allowed to remain wet for 5 minutes.
- When used on environmental, inanimate hard surfaces exhibits virucidal activity against Human Coronavirus and SARS associated Coronavirus when the treated surface is allowed to remain wet for 2 minutes.
- When used on environmental, inanimate hard surfaces exhibits virucidal activity against Rhinovirus type 39 (ATCC VR-340) and Rotavirus when the treated surface is allowed to remain wet for 3 minutes.
- When used on environmental, inanimate hard surfaces exhibits virucidal activity against HIV-1 (associated with AIDS) when the treated surface is allowed to remain wet for 1 minute (60 seconds).
- When used on environmental, inanimate hard surfaces exhibits effective virucidal activity against Norovirus (Norwalk virus), Feline Calicivirus and Rabies Virus when the treated surface is allowed to remain wet for 30 seconds.
- Kills HIV-1, HBV and HCV on pre-cleaned surfaces/objects

Read and understand the product's label and Material Safety Data Sheet ("MSDS") for precautionary and first aid information.

The MSDS is available on the Nu-Calgon website at [www.nucalgon.com](http://www.nucalgon.com) or is returnable by U.S. Mail upon request.

